Italiano

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EU regulatory agent:

US and international patents apply

Specifications are subject to change without notice

Sensors that make sense

User Manual

Respiration Pressure Flow Sensors

Table with 2 columns: Model, Connector/System. Rows include 1200, 1205, 1250, 1257, 1258, 1261, 1262, 1263, 1268, 1269.

Manufactured by: S.L.P. Ltd.
62 Anilevitz St, Tel Aviv, 67060 ISRAEL
Tel: +972 3 5371281, Fax: +972 3 5371282
www.sleepsense.com

Distributed by: S.L.P. Inc.
1121 E. Main St. #240, St. Charles, IL 60174, USA
Tel: (630) 513 7479, Fax: (630) 513 8478
E-mail: info@sleepsense.com
Toll-free(U.S.): 888-slp-sens (757 7367)

English

Thank you for purchasing a SleepSense sensor

Description

Snoring sensors respond to snoring and other sounds in the audio range picked up through the skin, and convert them to a small analog voltage that provides a clear, reliable indication of the presence of these sounds. Piezo or dynamic sensors do not require a battery or interface. Electret microphones are powered by the system. They are easily activated by plugging the output connector into the appropriate input on the head box.


- Dynamic snore sensors offer a wide frequency range, and sensitivity to other respiratory sounds besides snoring.
• Piezo sensors are very robust, and are relatively insensitive to ambient noises.
• Electret microphones offer very wide frequency response and very high sensitivity.

Indications for use

SleepSense® sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

Applying the Sensor

- The sensor is attached to the patient's neck.
• To position the sensor, place the sensor with the opening side or the protruding side (with the Piezo Sensor), against the skin and use medical grade tape (intended for facial skin) to hold the sensor in place. Make sure the sensor is tight against the skin to ensure that the vibrations are clearly transmitted to the sensor during snoring.
• Ask the patient to simulate snoring, and check that clear, strong signals are being transmitted before leaving the patient and starting the study.

Cleaning the Sensor

- Wipe the sensor and cable with a non-corrosive (to plastic) cleanser to clean before use.
• Make sure the complete sensor assembly is thoroughly dry before reusing it.
• To sterilize the sensors, use standard gas sterilization procedures.
• Do not autoclave or soak the sensor in disinfectants.

Technical Specifications:

For a complete list of technical specifications please contact us or visit our website.

Table with 4 columns: Sensor, Dynamic, Piezo, Electret. Rows include Cable, Sensor design, Recommended filter settings, Signal Output, Operating conditions, Storage temperature, and Operating and storage humidity.

Warnings and Precautions

- SleepSense® sensors are for professional use only.
• SleepSense® sensors are for sale by or on the order of a physician only.
• SleepSense® sensors may be used only in conjunction with an approved recording system.
• SleepSense® sensors are only intended for use on healthy skin.
• Discontinue use if sensor shows signs of wear, damage, or exposed metal.

Limited Warranty
SLP warrants the snoring sensor to be free of defects in materials and workmanship for a period of six months from the date purchased. The sole liability of SLP and our Distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period.

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