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Item N14XX - rev.2.0

Distributed by:

EU regulatory agent:  
Medes Limited  
9 Beaumont Hill, Radlett  
Herts WD7 7AR  
London, England  
Tel/Fax: +44 19238 59810

US and international patents apply

Specifications are subject to change without notice

**SleepSense®**

Sensors that make sense

**User Manual**

Thermal Flow Sensors

Model	Connector / System
1401	Thermocouple Flow sensor / Safety DIN connectors
1402	Thermocouple Flow sensor / 3.5mm
1405	Thermistor Flow sensor (Compumedics P Series™ )
1450B	Thermistor Flow sensor Interface box/ Safety DIN connectors
1450S	Thermistor Flow sensor
1468	Thermistor Flow sensor (short cable) / 2 pin safety connector
1469	Thermistor Flow sensor / 2 pin safety connector

**Manufactured by: S.L.P. Ltd.**  
62 Anilevitz St, Tel Aviv, 67060 ISRAEL  
Tel: +972 3 5371281, Fax: +972 3 5371282  
[www.sleepsense.com](http://www.sleepsense.com)

**Distributed by: S.L.P. Inc.**  
1121 E. Main St. #240, St. Charles, IL 60174, USA  
Tel: (630) 513 7479, Fax: (630) 513 8478  
E-mail: [info@sleepsense.com](mailto:info@sleepsense.com)  
Toll-free(U.S.): 888-slp-sens (757 7367)

# Reusable Flow Sensors

English

Thank you for purchasing a SleepSense sensor

## Description

Comprising three or four heat sensitive elements, the respiration flow sensor is placed under the patient's nostrils and produces a signal that is directly proportional to the temperature changes of air inhaled and exhaled during respiration. The sensor generates a small analog electrical voltage that provides a clear, reliable indication of respiration airflow. The sensors should be connected to the thermal-flow or to any AC (low level) channel on the recorder.

## Indications for use

SleepSense® sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

## Applying the Sensor

The respiration flow sensor is placed on the patient's upper lip, under the nostrils.

**The sensors are delicate and must be handled with care.**


## Positioning the Sensor

1. Bend the two nasal prongs forward to prevent them from entering the nostrils. Make a round bend slightly away from the sensor's body.
2. Align the tip of oral prong with the center of the mouth/upper lip.
3. Draw the wires on either side of the sensor over the patient's ears and use the clasp under the chin to secure the wires in place.

## Recommendations:

- Tape the wire on both cheeks to reduce strain and enhance reliability.
- Check that clear, strong signals are being transmitted before leaving the patient.

## Cleaning the Sensor

- Wipe the sensor and cable with a non-corrosive (to plastic) cleanser to clean before use.
- Make sure the complete sensor assembly is thoroughly dry before reusing it.
- To sterilize the sensors, use standard gas sterilization procedures.
- Do not autoclave or soak the sensor in disinfectants.

## Technical Specifications:

For a complete list of technical specifications please contact us or visit our website.

<b>Sensor:</b>	Thermocouple	Thermistor
<b>Technology:</b>	Small rounded TPE center with 3 prongs	
<b>Recommended filter settings:</b>	High pass: 0.1 Hz / Low pass: 70 Hz or as recommended in the system's instructions	
<b>Signal output:</b>	Approximately 0.2 mV for normal respiration	Approximately 1 mV for normal respiration with interface box
<b>Operating conditions:</b>	5°C (40°F) - 40°C (104°F)	
<b>Storage temperature:</b>	-20°C (-4°F) - 60°C (140°F)	
<b>Operating and storage humidity:</b>	5% - 95% (Non-condensing)	

## Warnings and Precautions

- SleepSense® sensors are for professional use only.
- SleepSense® sensors are for sale by or on the order of a physician only.
- SleepSense® sensors may be used only in conjunction with an approved recording system.
- SleepSense® sensors are only intended for use on healthy skin.
- Discontinue use if sensor shows signs of wear, damage, or exposed metal.

### Limited Warranty

SLP warrants the respiration-flow sensor to be free of defects in materials and workmanship for a period of six months from the date purchased. The sole liability of SLP and our distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period. Under no circumstances shall SLP or our distributor(s) be liable for any loss of revenues or damage, direct, consequential, or incidental, including loss of profit, property damage, or personal injury arising from the use of, or the inability to use this product. This warranty is intended only for the original buyer and is in lieu of all other warranties or previous agreements, expressed or implied. This warranty is rendered void if the product is used for other than its intended purpose or is subject to abuse, misuse, tampering, neglect, or unauthorized modifications. Use of this product constitutes acceptance of this warranty in total.

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