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Better sensors. Better sense.

User Manual

Inductive Plethysmography Effort Sensors (White/Blue)

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www.sleepsense.com

English

Thank you for purchasing a SleepSense® sensor.

### Description

The inductance of the inductive bands is exactly proportional to the cross section area of the band's closed loop around the body. It therefore changes as the patient breathes. The electronics in the system or interface box measure these inductance changes to generate a signal that is the most accurate and reliable replication of respiration waveforms. Inductive bands easily connect to any system compatible with inductive effort sensors through a dedicated interface cable. The SleepSense Inductive Interface Cable allows the user to enjoy the benefits of inductive effort sensing with other systems as well.

### Indications for use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

### Applying the Sensor

The patient should lie on his/her back for optimal positioning of the sensor.


### Positioning the Sensor

- Select the proper size band for the patient.
- 1. Strap the sensor bands around the patient's abdomen or chest, over any clothing.
- 2. Attach the first wire to the snap that is on the same side of the band as the Velcro® patch.
- 3. While holding this end to the left of the patient's centerline, wrap the band around the patient's chest or abdomen. For smaller patients move the starting point further to the left.
- 4. Pull the other end over the first end and secure to the band on the right side of the centerline. The band should be snug but not too tight, (see image).
- 5. Attach the second wire to the exposed snap. Secure both wires about 3" away from the connection.
- Check that clear, strong signals are being transmitted before leaving the patient.
- **When removing the sensor band do not pull on the wires. Remove from the snap.**

Attention! Using two chest or two abdomen interface cables on the same patient may cause interference.

### Cleaning the Sensor

- The bands are machine washable.
- Wipe the Interface Cable with a non-corrosive (to plastic) cleanser to clean before use.
- Make sure the bands and Interface Cable are thoroughly dry before reusing it.
- To sterilize the band / Interface Cable, use standard gas sterilization procedures.
- Do not autoclave or soak the Interface Cable in disinfectants.

### Technical Specifications:

For a complete list of technical specifications please contact us or visit our website.

Item #:	9003/9002 (different lengths)	Interface Cables
<b>Sensor:</b>	Synthetic inductive bands - non latex	Inductive Effort Sensor interface cable
<b>Recommended filter settings:</b>	High pass: 0.1 Hz / Low pass: 40 Hz or as recommended in the system's instructions	
<b>Operating conditions:</b>	5°C (40°F) - 40°C (104°F)	
<b>Storage temperature:</b>	-20°C (-4°F) - 60°C (140°F)	
<b>Operating and storage humidity:</b>	5% - 95% (Non-condensing)	

### Warnings and Precautions

- SleepSense sensors are for professional use only.
- Caution: US Federal Law restricts this device to sale, distribution and use by or on the order of a physician.
- SleepSense sensors may be used only in conjunction with an approved recording system.
- SleepSense sensors are only intended for use on healthy skin.
- Discontinue use if sensor shows signs of wear, damage, or exposed metal.

**Attention:** contains a lithium battery - dispose of properly.

Limited Warranty  
SLP warrants the Inductive Interface cable to be free of defects in materials and workmanship for a period of twelve months from the date purchased. The sole liability of SLP and our Distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period. Under no circumstances shall SLP or our Distributor(s) be liable for any loss of revenues or damage, direct, consequential, or incidental, including loss of profit, property damage, or personal injury arising from the use of, or the inability to use this product. This Warranty is intended only for the original buyer and is in lieu of all other warranties or previous agreements, expressed or implied. This warranty is rendered void if the product is used for other than its intended purpose or is subject to abuse, misuse, tampering, neglect, or unauthorized modifications. Use of this product constitutes acceptance of this warranty in total.

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