Thank you for purchasing a SleepSense® sensor.

Description

The inductance of the inductive bands is proportional to the cross section area of the band's closed loop around the body. It therefore changes as the patient breathes. The electronics in the interface box measure the thoracic and abdominal effort signals and their sum (sum model) to provide accurate and reliable respiratory effort signals used to assist the technician in identifying desynchronized chest and abdominal breathing.

Instructions for Use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

Applying the Sensor

The patient should lay on his/her back for optimal positioning of both bands.

For use with SLP inductive bands.

Positioning the Sensor (See Figure A)

- 1. Strap the appropriately sized sensor bands around the patient's chest and abdomen, over any clothing.
2. Attach the first wire to the snap that is on the same side of the band as the Velcro® patch.
3. While holding this end to the left of the patient's centerline, wrap the bands around the patient's chest and abdomen.
4. Pull the other end over the first end and secure to the band on the right side of the centerline.
5. Attach the second wire to the exposed snap.
When removing the sensor band do not pull on the wires. Remove from the snap.

Interface Box Setup (See Figure B)

- The chest interface cable MUST be connected to chest inputs to activate the unit.
Once connected, a green light appears indicating the device is ON.
A red light appears when the interface box is operating and the battery power is running low.
When this occurs, it is recommended to replace the batteries within the next few days.
Plug the black connector of the chest interface cable into the designated socket on the interface box.
Plug the black connector of the abdomen interface cable into the designated socket on the interface box.
Plug the black output connector into the designated socket on the interface box.
When removing the interface cables do not pull on the wires. Remove from the connector.
Check that clear, strong signals are being transmitted before leaving the patient.

Cleaning the Sensor

- The bands are machine washable.
Wipe the Interface Cable with a non-corrosive (to plastic) cleanser to clean before use.
Make sure the bands and Interface Cable are thoroughly dry before reusing it.
To sterilize the band / Interface Cable, use standard gas sterilization procedures.
Do not autoclave or soak the Interface Cable in disinfectants.

Technical Specifications:

For a complete list of technical specifications please contact us or visit our website.

Table with 2 columns: Items #, and specifications for Box 9201, Chest cable PS9203, Abdomen Cable PS9202, Output cable PS9204. Includes filter settings, operating conditions, storage temperature, and humidity.

Warnings and Precautions section containing safety instructions for SleepSense sensors, including professional use only, caution about federal law, and attention to battery disposal.

Limited Warranty
SLP warrants the Multi-Rip box and cables to be free of defects in materials and workmanship for a period of twelve months from the date purchased.

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Table with 2 columns: Article #, and specifications for Boîte 9201, câble de poitrine PS9203, câble d'abdomen PS9202, Câble de sortie PS9204. Includes filter settings, operating conditions, storage temperature, and humidity.

Avertissements et Précautions section containing safety instructions in French, including professional use only, caution about federal law, and attention to battery disposal.

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Table with 2 columns: Artículos #, and specifications for Caja 9201, Cable de tórax PS9203, Cable de abdomen PS9202, Cable de salida PS9204. Includes filter settings, operating conditions, storage temperature, and humidity.

Advertencias y precauciones section containing safety instructions in Spanish, including professional use only, caution about federal law, and attention to battery disposal.

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